FUS-3000Plus NEW

Urinalysis Hybrid



Technical Speci	fication
Sample Type	Urine
Test Items	Formed Elements Items: NRBC, MIRBC, ARBC, SRBC, OSRBC, WBC, WBCC, SQEP, RTEP, TREP, HYAL,
	GRAN, WAXY, BROAD, OCAS, BACI, SUCO, HYST, BYST, CAOX, URIC, MAPH, OCRY, SPRM, MUCS
	Chemistry Items: UBG, BIL, KET, PRO, BLD, LEU, NIT, GLU, pH, SG, VC(Optional), MALB(Optional),
	CRE(Optional), CA(Optional).
	Physics Items: Conductivity(Optional), Specific Gravity(Optional), Color(Optional), Turbidity (Optional)
Principles	Formed Elements: Flow Imaging Technique
	Chemistry: Photoelectric Colorimetry
Sample loading	50/270 Samples (Optional)
Throughput	Hybrid Mode/ Formed Element Mode: 120 Tests/ Hour
	Chemistry Mode: 240 Tests / Hour
Min Sample Volume	Hybrid mode / Formed Elements Mode: 3.0 mL, Chemistry mode: 1.5 mL
Aspiration Volume	Hybrid mode / Formed Elements Mode: 2.2 mL, Chemistry mode: 0.4 mL
Interface	Communication: RS-232 Port; LIS: Serial Port & Net Port
LIS/ HIS	Bidirectional communication with LIS/HIS
Bar Code Identification	n Internal Bar Code Scanner
	External Bar Code Scanner(Optional)
Memory Capacity	≥ 100,000
Dimension	$779 \text{mm} \times 688 \text{mm} \times 584 \text{mm} (\text{L} \times \text{W} \times \text{H}), \ 897 \text{mm} \times 1566 \text{mm} \times 584 \text{mm} (\text{L} \times \text{W} \times \text{H}) \ \text{with pre-storage \& reclaim tray}$
Weight	82 kg
Operation Temp.	10-30°C
Humidity	≤75%
Power	100~240V 50/60Hz

DIRUI INDUSTRIAL CO.,LTD.

·Specifications subject to change without notice.





3333 Yiju Street,New&High Tech. Development Zone Changchun,Jilin 130103,P.R.China Tel:+86(431)81935331 85100409 Fax:+86(431)85172581 85083741 E-mail:dirui@dirui.com.cn Http://www.dirui.com.cn









THE LEADER IN URINALYSIS



FUS-3000Plus

■ Urinalysis Hybrid

- Inherit classical flow imaging analysis technology
- Integrate physical, chemical and formed element analysis
- FUS-3000Plus—The leader in Urinalysis



Advanced Principle

Comprehensive Parameters

Efficient Test
Stable System



ADVANCED PRINCIPLE

The FUS-3000Plus utilizes flow imaging analysis technology for accurate and reliable results

Sheath Flow Technology


The urine formed element is directed by the sheath liquid with a single layer to avoid gather or overlap.

High-Speed Photograph Technology

2000 frames image are taken with 400 times magnification to increase detection rate.

Artificial Intelligent Identification(All)

All makes the result highly concordant with Golden Standard.



Identification and classification





COMPREHENSIVE PARAMETERS

FUS-3000Plus incredible provides the most comprehensive 43 report parameters and 2 research parameters results in the world

Formed Element Parameter

- 25 formed element parameters breakthrough the technique difficulties of urine formed element automatic analysis.
- Particles sub-classification offer high valuable reference for clinic diagnosis

Abnormal RBC







Abnormal RBC alert uniformity, non-uniformity or miscibility RBC Abnormal RBC ratio alert the hematuresis resource

Epithelial Cell









Epithelial Cell indicate the lesion position of urinary system diseases, such as urinary tract inflammation, nephric tubule acute necrosis and urocystitis

Pathological CAST







Pathological CAST indicate kidney parenchymal damage such as glomerulonephritis, chronic kidney failure, etc

Crystal









Crystal sub-classification indicate kidney stone and may indicate metabolic disorders

Bacterial







Bacterial is classified and identified
Urinary Tract Infection(UTI) is alerted based on bacterial,
WBC, RBC and chemistry result


Chemical and Physical Parameter

- 14 chemical parameters + 4 physical parameters
- ACR offers a precise estimate of the 24 hours urinary albumin excretion and the advanced detection of early-stage renal dysfunction
- Conductivity indicate the concentration and dilution ability of Kidney





Complete quality assurance system Complete calibrator, quality control, strip & reagent Specific complex quality control Build up stable qualified test system





Sieve test function: intelligently activate formed element test



One key operation: execute all works with only one key